

LABORATORY TROLLEYS



AMADWALA LABORATORY TROLLEYS are ACID (CORROSIVE) resistant and impervious to Acids and most liquids. Manufactured entirely with high Acid / Corrosivity Resistant 10mm thick Polypropylene

FEATURES:

- Fully Acid/Corrosive Resistant Polypropylene Structure and Shelves
- Impervious to most liquids
- Fully Thermo Welded Construction – No leaks
- Dished Shelving – Spillage containment
- Two (2) or Three (3) Shelf units available
- Lockable Swivel Castors x 4
- Hard wearing

PECIFICATION

MODEL	CONFIGURATION	SIZE (mm)			SHELF QTY
		Height	Length	Width	
2LT	Two Shelf	945	770	600	2
3LT	Three Shelf	945	770	600	3

**** Special Sizes can be manufactured on request**

LITY

	(mm)	
MODEL	2LT	3LT
A	945	945
B	853	853
C	710	710
D	770	770
E	750	750
F	600	600
G	580	580
H	523	255

