

FIRECLAY CRUCIBLES NOMINAL VOLUMES AND DIMENSIONS



<u>FIRED SIZE</u>	<u>Height</u>	<u>Top D</u>	<u>Bottom D</u>	<u>Volume</u>	<u>Package</u>
No 3	110	90	64	300	1000
No 4	120	115	70	600	1000
No 5	155	115	74	650	800
40g	135	92	58	393	1550
40gT/45T	160	90	58	450	1500
50g	142	104	59	500	1500
65g	165	105	60	600	1200

PACKED IN CRATES

MANUFACTURED IN THE REPUBLIC OF SOUTH AFRICA