

MAGESIUM OXIDE BASED CUPELS NOMINAL WEIGHTS AND DIMENSIONS



SIZE
S.

<u>CUPELS</u>	<u>WIDTH</u>	<u>HEIGHT</u>	<u>DEPTH</u>	<u>WEIGHT</u>	<u>PACK</u>
BULLION					
6A	31.8	28.2	12	50G	300
7A	31.8	28.2	9.5	60G	300
8	44.6	30.2	11.8	80G	500
9	50.3	30.2	14.3	98G	500
10	60.0	47	8.1	183G	250

- Features of these CUPELS are:
- Extensively tested : Known for their Hardness and Durability :
- No lead spitting :
- Rapid absorption of lead.
- These CUPELS are manufactured with high grade Magnesium oxide.
- Low losses.
- Free from cracking.

**ALL CUPELS ARE PACKED 48 BOXES PER PALLET 20 PALLETS
PER CONTAINER:**

MANUFACTURED IN THE REPUBLIC OF SOUTH AFRICA

CUPELS	PER BOX	WIDTH	HEIGHT	DEPTH	WEIGHT
BULLION 500		25.4	26.6	7.1	27GRM
6A 300		31.8	28.2	12.0	50GRM
7A 300		31.8	28.2	9.5	60GRM
7R 500		34.9	34.9	9.0	53GRM
8 500		44.6	30.2	11.8	80GRM
9 500		50.3	30.2	14.3	98GRM
9S 500		50.8	38.2	14.3	122GRM
10 250		60.0	47	8.1	183GRM
BLOCK 30 500		54	54	24 (block)	

MUFFLE DRESSING 50KG BAG FINE

REVERBERITE 50KG BAG

Reg No: 2013/084255/07
Vat No: 4410110631
Telephone: 011-434-1450
Unit D1
Selby Industrial Park,
Cnr John and Prop Street,
Selby, Johannesburg

**YOUR PARTNER IN ADVANCED
MINING LABORATORY SOLUTIONS**

